F

**KELO-COTE®  
UV**

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**KELO-COTE®  
UV**

Advanced formula scar gel

**Kelo-cote® UV** with SPF 30 was specifically developed to protect new and old scars from the aggravating effects of the sun by preventing sun-induced redness and hyper-pigmentation while improving overall appearance. It can be used on all scars resulting from surgery, trauma, wounds and burns. **Kelo-cote® UV** patented silicone technology has been shown in clinical studies to flatten, soften and smooth scars, as well as relieve associated redness, itch and pain. **Kelo-cote® UV** may also help in the prevention of excessive and abnormal scar formations. This unique technology allows the gel to self-dry and form a flexible, breathable, waterproof layer that binds to the stratum corneum (the outer layer of skin), providing hydration and protection to the scar site.

Whether a scar is a result of a sutured incision or the natural healing of a wound, it should be protected from the sun during its whole maturation process. (The time that it remains red or pink), which can take from 6 months to 3 years according to studied cases. If a fresh or recent scar is left unprotected from the sun, it may result in unwanted long-term pigmentation of the scar. This increased pigmentation can be aesthetically compared to the surrounding normal skin. Therefore, it is recommended that a protective sunscreen be applied to scars, and be reapplied every 2 or 3 hours when the body is exposed to the sun.\*\*

**Kelo-cote® UV**'s high protection sunscreens and patented, proven silicone technology ensures the optimal healing environment for all scars. This unique blend of protective and treatment ingredients in one tube offers the convenience of a single application and promotes an aesthetic scar outcome.

**Kelo-cote® UV** is ideal for the face, neck, hands, chest or any other area exposed to the sun. **Kelo-cote® UV** can be used on joints, flexures, hairy areas and large area scars such as those caused by burns or major trauma. **Kelo-cote® UV** may also be used in conjunction with pressure garments or other scar treatments.

#### Intended use:

**Kelo-cote® UV** is a unique scar therapy providing a patented silicone gel to help prevent abnormal scar formation with high protection sunscreen agents. For maximum effect, **Kelo-cote® UV** should have 24 hour contact with the skin. Recommended minimum treatment time is 60-90 days. Larger and older scars may take longer and continued use is recommended if improvement is still seen after the initial 90 days.

Ingredients	Purpose
Octinoxate 7.5%	.....Sunscreen
Octisalate 5.0%	.....Sunscreen
Octocrylene 10%	.....Sunscreen
Oxybenzone 6.0%	.....Sunscreen
Polysiloxanes	
Silicone Dioxide	
Isopropyl Myristate	

Contents: 6g & 15g

#### Warnings:

For external use only

• Avoid direct contact with eyes, mucous membranes, third degree burns and open wounds.

• Should not be used on dermatological conditions that disrupt the integrity of the skin.

• If irritation occurs, discontinue use and consult your physician.

• Keep out of reach of children.

#### Directions:

• Ensure that the affected area is clean and dry.

• Apply evenly in a thin layer twice daily, 30 minutes before sun exposure.

• Once dry, cosmetics may be applied.

• Keep babies and young children out of direct sunlight.

#### Tips for drying:

If **Kelo-cote® UV** has not dried within 4-5 minutes, you have probably used too much. Gently remove the excess and allow the drying process to continue.

**Kelo-cote® UV**'s drying process is assisted by the temperature difference between the gel and the skin. To reduce drying time in hotter climates, keep **Kelo-cote® UV** in the refrigerator. In colder weather, use the low setting on a hair dryer to reduce drying time.

#### Other Information

• Store below 25°/77°F

www.kelocote.co.uk

Device Authorised Representative Alliance Pharma S.r.l. Viale Restelli 5, 20124 Milano, Italy

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☒ Bevollmächtigter Vertreter Alliance Pharma S.r.l. Viale Restelli 5, 20124 Milano, Italy
Advanced Bio-Technologies, Inc 4830 West Kennedy Blvd, Suite 600, Tampa, Florida 33609 USA